

Product Information

pro.Tray – Ultra fine baskets

Excellent solution for small parts



pro.Tray – Ultra fine baskets, Art.-Nr.: PT-034-A2

pro.Tray

Ultra fine baskets are used for small parts, delicate micro surgical and laboratory instruments for cleaning, disinfection and sterilization.

- ✓ Excellent protection
- ✓ Excellent cleaning
- ✓ Excellent drying

Innovative solution for processing of medical devices.
Patient safety and protection for your investment.

- please note back page -

pro.Tray – Ultra fine baskets

Ultra fine mesh with opening square 0.5 mm and wire diameter 0.3 mm.
Top and lid frame wire diameter 3.0 mm.
Wire locking clip.

Size in mm	Art.-Nr.
80 x 40 x 30	PT-034-A1
80 x 80 x 34	PT-034-A2
100 x 100 x 50	PT-034-B
180 x 130 x 20	PT-034-C
230 x 100 x 60	PT-034-D
240 x 120 x 60	PT-034-E

For further information, questions or orders, please visit our website [www.med**proces**stec.com](http://www.medprocesstec.com) or contact us:

[contact@med**proces**stec.com](mailto:contact@medprocesstec.com)